

CE
0473

PORTEX™

en acapella® Vibratory PEP Therapy System

fr acapella® Système de thérapie par PEP vibratoire

de acapella® PEP-Vibrationstherapiesystem

it Dispositivo PEP a vibrazioni acapella®

es Sistema acapella® de PEP vibratoria

pt Sistema de Terapia de Pressão Expiratória Positiva (PEP) Vibratório acapella®

nl acapella® Vibratie-PEP-therapiesysteem

sv acapella® vibrerande system för PEP-behandling

da acapella® PEP-system med vibration.

Smiths Medical ASD, Inc.
Keene, NH 03431, USA

Smiths Medical International Ltd.
Hythe, Kent, CT21 6JL, UK

Australian Representative:
Smiths Medical Australasia Pty. Ltd.
Brisbane, QLD 4113, Australia

www.smiths-medical.com


600-442-004 REV.002

smiths medical



Figure 1.
Figure 1.
Abb. 1.
Figura 1.
Figura 1.

Figura 1.
Figuur 1.
Figur 1.
Figur 1.



<p>Counterclockwise Decreases frequency-resistance.</p> <p>Dans le sens inverse des aiguilles d'une montre - Diminue la fréquence-résistance.</p> <p>Drehung entgegen dem Uhrzeigersinn verringert den Frequenzwiderstand.</p> <p>In senso antiorario per ridurre resistenza-frequenza.</p> <p>De dcha. a izda. reduce la resistencia-frecuencia.</p> <p>A rotação no sentido contrário aos ponteiros do relógio diminui a frequência-resistência.</p> <p>Naar links draaien verlaagt de frequentieweerstand.</p> <p>Motus minskar frekvensen/motståndet.</p> <p>Mod uret mindsker frekvensmodstand.</p>	<p>Clockwise Increases frequency-resistance.</p> <p>Dans le sens des aiguilles d'une montre - Augmente la fréquence-résistance.</p> <p>Drehung im Uhrzeigersinn erhöht den Frequenzwiderstand.</p> <p>In senso orario per aumentare resistenza-frequenza.</p> <p>De izda. a dcha. aumenta la resistencia-frecuencia.</p> <p>A rotação ao sentido dos ponteiros do relógio aumenta a frequência-resistência.</p> <p>Naar rechts draaien verhoogt de frequentieweerstand.</p> <p>Med ur et ökar frekvensen/motståndet.</p> <p>Med uret öger frekvensmodstand.</p>
--	---

Figure 2.
Figure 2.
Abb. 2.
Figura 2.
Figura 2.

Figura 2.
Figuur 2.
Figur 2.
Figur 2.

en English

These instructions contain important information for safe use of the product. Read the entire contents of these instructions for use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient. It is the responsibility of the healthcare practitioner to assure that the instructions for use and maintenance are understood by and provided to the care giver.

1. **DESCRIPTION:**
The acapella® vibratory PEP system is a single patient use device that provides Positive Expiratory Pressure (PEP) therapy for patients who have Cystic Fibrosis, COPD, asthma, and lung diseases with secretory problems, and patients with atelectasis. All patients must be capable of following instructions for Positive Expiratory Pressure Therapy. Review the diagram and the product to become familiar with all of the product features.
The acapella® system consists of:
1. Removable mouthpiece (mask options are also available)
2. acapella® vibratory PEP device with 22 mm male fittings (See Figure 1) including:
A. Expiratory resistance/frequency adjustment dial
B. One-way inspiratory valve
C. Patient end which adapts to a mouthpiece, or mask. (If pressure monitoring feedback is desired, the Smiths Medical TheraPEP® pressure port may be inserted (See Figure 4).
D. Mouthpiece (detachable)
3. **INDICATIONS:**
The Smiths Medical acapella® is intended for use as a Positive Expiratory Pressure (PEP) device. It may also be used simultaneously with nebulized aerosol drug delivery.
4. **CONTRAINDICATIONS:**
Although no absolute contraindications to the use of PEP Therapy have been reported, the following should be carefully evaluated before a decision is made to initiate therapy:
• Inability to tolerate increased work of breathing
• Hemodynamic instability
• Intracranial pressure (ICP) > 20 mm Hg
• Acute sinusitis
• Recent facial, oral or skull surgery or trauma
• Epistaxis
• Esophageal surgery
• Active hemoptysis
• Untreated pneumothorax
• Nausea
• Known or suspected tympanic membrane rupture or other middle ear pathology
5. **WARNINGS:**
Use of this device at excessive pressures may have adverse effects. Expiratory pressures above 20 cm H₂O in patients sensitive to increased transpulmonary pressure may develop one or more of the adverse side effects listed below.
Expert clinical judgment should be exercised in the selection of the appropriate setting for each individual patient. Failure to match the appropriate resistance setting on the +/- dial indicator with the patient's expiratory flow may result in failure to achieve therapeutic objectives of vibratory PEP therapy or one or more adverse side effects below.
Adverse reactions may include:
• Increased work of breathing that may lead to hypoventilation and hypercarbia
• Increased cranial pressure
• Cardiovascular compromise
• Myocardial ischemia
• Decreased venous return
• Air swallowing with increased likelihood of vomiting and aspiration
• Claustrophobia
• Skin break down and discomfort from mask
• Pulmonary barotraumas
6. **PRECAUTIONS:**
6.1 Bleach is not recommended for use on the acapella®. It may deteriorate the nickel plated mechanism located in the interior of the device.
6.2 DO NOT MICROWAVE. The metal and magnet might ignite.
6.3 It is the responsibility of the user to ensure all sterility verification(s).
6.4 Visually inspect the device to ensure that the unit is free of contamination and foreign objects.
6.5 Verify all connections are secure.
7. **SUGGESTED INSTRUCTIONS FOR USE**
7.1 **Initial Settings:**
6.1.1 If this is the first use of the acapella® device, ensure that the frequency adjustment dial is turned counterclockwise to its lowest frequency-resistance setting. (See Figure 2)
6.1.2 Instruct the patient to relax while performing diaphragmatic breathing and inspiring a volume of air larger than normal tidal volume (but not to total lung capacity).
6.1.3 Direct the patient to exhale to Functional Residual Capacity (FRC) actively, but not forcefully, through the device.
6.1.4 The patient should be able to exhale for 3-4 seconds while the device vibrates. If the patient cannot maintain an exhalation for this length of time, adjust the dial clockwise (See Figure 2). Clockwise adjustment increases the resistance of the vibrating orifice, which will allow the patient to exhale at a lower flow rate.
6.1.5 Selection of the proper resistance range produces the desired inspiratory to expiratory (I:E) ratio of 1:3 to 1:4.
Note:
If the desired resistance and I:E ratio cannot be achieved, consider using an acapella® device designed for alternate flow ranges.
6.1.6 Once the proper range has been identified, the patient may be instructed to exhale harder or softer, or dial adjustments may be made to optimize the response the user "feels" from the vibratory pressure. Several uses may be needed to ensure that individual patient needs are being met.
6.1.7 **Procedure for the User:**
6.2.1 Ensure the adjustment dial is set to the correct range as identified by your clinician.
6.2.2 Sit with elbows resting comfortably on table.
6.2.3 Place mouthpiece lightly in mouth.
• Be sure to maintain a tight seal on the mouthpiece during exhalation.
• Your clinician may recommend the use of a nose clip, if necessary.
• If using a mask, apply mask tightly but comfortably over nose and mouth.

- 6.2.4 Breathe from the diaphragm, as directed by your clinician, taking in a larger than normal breath, but not filling your lungs to capacity.
- 6.2.5 Hold your breath for 2-3 seconds.
- 6.2.6 Exhale actively, but not forcefully, through the device. Exhalation should last approximately 3 to 4 times longer than inhalation.
- 6.2.7 Perform 10-20 PEP breaths as recommended by clinician.
- 6.2.8 Remove mouthpiece (mask) and perform 2-3 "huff" coughs to raise secretions as needed. Your clinician may direct you on proper cough technique.
- 6.2.9 Repeat steps 6.2.2 to 6.2.7 as prescribed.
Note: See nebulizer setup section.
- 6.3 **Set-up Nebulizer**
6.3.1 Review the diagrams contained with the device.
6.3.2 A possible nebulizer and acapella® setup is reflected below. (See Figure 3)
6.3.3 Follow setup instructions for each device.
6.3.4 Follow cleaning instructions contained with each device.
6.3.5 Inspect device(s) on a routine basis to ensure proper use and function.
6.3.6 If damaged, do not use.
6.3.7 Verify all connections are secure.
- 6.4 **Set up pressure gauge/indicator** (See Figure 4)
6.4.1 Place the pressure gauge/indicator (Smiths 20-0010) between the mouthpiece and the device.
6.4.2 Verify all connections are secure.
7. **CLEANING AND DISINFECTING INSTRUCTIONS:**
Precaution: Bleach is not recommended for use on the acapella®. It may deteriorate the nickel plated mechanism located in the interior of the device.
7.1 **Cleaning: This should be done prior to Disinfecting (7.2)**
As per the Cystic Fibrosis Foundation's cleaning and disinfecting guidelines entitled, "Respiratory: Stopping the Spread of Germs" 2008, below are the guidelines for acapella®.
Cleaning with Liquid Dish Detergent:
As needed, detach mouthpiece (mask) then soak the device and mouthpiece in warm, soapy water as required to remove visible contaminants. Use a liquid dish detergent (Dawn or equivalent), mixing two (2) tablespoons of detergent per one (1) gallon water. Rinse thoroughly with sterile water, and allow parts to air dry. Drain the device by placing it in a normal resting position. (See Figure 1)
7.2 **Disinfecting:**
8. **DISPOSAL:**
Dispose of acapella® in a safe manner according to Federal/State/Local regulations and guidelines for disposal of contaminated medical waste.
9. **ESTIMATED DEVICE LIFETIME:**
Assuming the manufacturer recommended cleaning protocol is followed, the acapella® device should have a useful life of six (6) months under normal and customary usage. The six (6) month usage duration is measured from the date of initial use.

REF	acapella® DM (Blue) Catalog Number	acapella® DH (Green) Catalog Number	Quantity Per Case
Mouthpiece	21-1015	21-1330	10
Pediatric Mask	21-3015	21-3330	1
Medium Mask	21-5015	21-5330	1
Large Mask	21-7015	21-7330	1

Accessory Products	Catalog Number	Catalog Number	Quantity Per Case
Pressure Port	20-0010	20-0010	10
TheraPEP pressure indicator, tubing, and pressure port	20-0022	20-0022	10

The Smiths Medical and Portex design marks and acapella are trademarks of the Smiths Medical family of companies. The symbol "i" indicates the trademark is registered in the U.S. Patent & Trademark Office & certain other countries. The products described are covered by one or more of the following U.S. Patent Nos. US 7,059,324, US 6,381,998 and US 440,651; other patent(s) pending. © 2010 Smiths Medical family of companies. All rights reserved.




Figure 3.
Figure 3.
Abb. 3.
Figura 3.
Figura 3.

Figura 3.
Figuur 3.
Figur 3.
Figur 3.




Figure 4.
Figure 4.
Abb. 4.
Figura 4.
Figura 4.

Figura 4.
Figuur 4.
Figur 4.
Figur 4.

acapella® with TheraPEP pressure indicator, tubing and pressure port.
acapella® avec indicateur de pression, tubulure et port de pression TheraPEP.
acapella® mit TheraPEP-Druckanzeige, Schläuchen und Druckanschluss.
acapella® con indicatore di pressione TheraPEP, tubo e porta di pressione.
acapella® con indicador de presión TheraPEP, tubo y puerto de presión.
acapella® com indicador de pressão TheraPEP, tubo e porta de pressão.
acapella® met TheraPEP drukindicator, leiding en drukpoort.
acapella® med TheraPEP tryckindikator, slang och tryckport.
acapella® med TheraPEP trykmåler, slange og trykport.